

SleepAngel®



SleepAngel® Express Pillows & Duvets

Cleaning and Care Instructions

SleepAngel® devices can be put into service immediately.
The product is not sterilized

1. Inspection and maintenance

Product inspections

- Visually examine the products to identify any breach of cover or significant damage.
- Replace any damaged products.
- When using wipes on all polyurethane surfaces, the surface should be wiped off using clean water after the wipe has been applied.

Weekly: (or on patient discharge)

- Check for signs of fluid contamination on the product and damage or holes in the product cover.
- Dispose of contaminated product as hospital waste.

2. Cleaning Instructions

For normal cleaning: use mild detergent and water solution. If the cover is heavily soiled, a 0.1% Chlorine solution (1,000 ppm) and cold water can be used. If the cover has been exposed to bodily fluids such as blood, a 1% Chlorine solution (10,000 ppm) and cold water can be used. Frequent or prolonged use of a 1% Chlorine solution (10,000 ppm) may reduce the useful life of fabric.

Note. Use of alcohol wipes and cleaning solutions may make a temporary change to the surface of the cover fabric and leave more susceptible to mechanical damage. Allow to dry fully before reuse. After processing, replace sheet carefully to help prevent mechanical damage.

Same cleaning procedure is valid for the PneumaPure™ filter. Avoid abrasive scrubbing of the filter and fabric.



Do not bleach, do not machine wash, do not dry clean, do not tumble dry, do not iron.

3. Storage

- The pillow can be stored flat, on the side, or on end without damage. Do not bend the pillow.
- The duvet can be stored flat folded or rolled, avoid bending and stretching the PneumaPure™ filter.



SleepAngel® Express Pillows & Duvets

Cleaning and Care Instructions

SleepAngel® devices can be put into service immediately.
The product is not sterilized

Cleaning / disinfectant products that can be used:

- % 1 Hydrogene Peroxide Solution
- Amonium Compounds
- Anios Oxy'Floor
- Aniospray Quick
- Anioxy Spray WS
- AntiGone
- Bacillol AF
- Bacillol Foam
- Bactacid Wipes
- Cavi Cide Surface Decontaminant Cleaner
- Cavi Wipes
- Chemisept wipes
- Chlorine Dioxide Solution, 1000 ppm
- Chlorine Releasing Solution, 1000 ppm
- Clinell Alcohol Wipes
- Clinell Sporicidal Wipe
- Clinell Universal Spray
- Descosept PUR Wipes RTU
- Desinfektionstuch
- Ecopine, (%5 Solution in Water)
- Ecopine, (%5 Solution in Water)
- Incidin Foam & Incidin Plus
- Incidin Pro
- Isopropyl Alcohol, %70
- Kohrsolin FF
- Life Clean
- Meliseptol Foam Pure
- Meliseptol Rapid
- Microbac Forte
- Mikrozyd AF Liquid

- Mikrozyd Sensitive Liquid
- NaOH Soluion
- Oxivir Plus
- PDI Sani Cloth AF3 Wipes
- PDI Sani-Cloth Bleach Germicidal Disposable Wipe
- Perform 1% and 3%
- Phenol Solution, %1
- Reinigungstuch cleaning wipes
- Rely+ON Perasafe
- Rely+ON Virkon
- Sani-Cloth Active
- Spray Off
- Sterisept wipes
- Super Sani Cloth Germicidal Disposable Wipe
- Terralin® protect (max. 4%)
- TPH protect
- Virusolve (%5 & %10 Solution in Water)
- Wet Wipe

Cleaning products that are not recommended and can hurt fabric surface:

- Chlorine Dioxide Solution, 10.000 ppm
- Disclean ENDO PAA (Peracetic Acid Solution)
- PDI Sani-Cloth Bleach Germicidal Disposable Wipe
- Peracetic Acid
- Incidin OxyFoam S

This list is not exhaustive and the appropriate concentrations should be strictly followed and not exceeded